

PACS Settings

30/07/2025 4:27 pm BST

Summary

PACS Settings can be accessed via: Administration > Settings and selecting the 'PACS' category.

This page documents each setting.

Details

The screenshot shows the 'Configuration Settings' interface. On the left, a vertical list of categories is displayed, with 'PACS' highlighted in a red box. The main area on the right is titled 'Picture Archiving and Communication System' and contains a dropdown menu. At the top, the 'Site' is set to 'Esher SA' and there is a 'Save' button.

Categories	Site	Save
Account	Esher SA	Save
Reporting		
Diary		
Reception		
Treatment		
Referral		
Lab		
Vet Scan Lab		
Labelling		
SMS		
Email		
DayLists		
System		
Miscellaneous		
HTML Labelling		
Smart Flow		
VetCheck		
Follow Up		
GDPR		
PACS		
Message		

Picture Archiving and Communication System

[Dropdown Menu]

Picture Archiving and Communication System

Specifies the PACS system used for the site selected. The list is produced from the PACS accounts added into Administration > PACS.